



PANAMEDICAL



# LINHA DE PRODUTOS

DÚVIDAS E ORÇAMENTOS

[comercial@panamedical.com.br](mailto:comercial@panamedical.com.br) | 11 23441900



CP ANGLE DISSECTOR



**CP-6 TUMOR DISSECTOR**  
 TL:19.5 cm  
 WL: 8.5 cm  
 Tip Size: 1 mm  
 Sharp  
 Round  
 07-821-06



**CP-7 TUMOR DISSECTOR**  
 TL:19.5 cm  
 WL: 8.5 cm  
 Tip Size: 1 mm  
 Curved  
 Round  
 07-821-07

*Shown above at 100% scale  
 Details shown larger*

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

Part numbers may change without notice.  
 For US commercial use only.  
 Refer to the Instruction for Use before using.  
 The colors in the photos may differ from the actual items.

[www.mizuho.com](http://www.mizuho.com)



©Mizuho America, ABPM 005, Rev A, Effective Date 09/28/2018

SI-3-1

DISSECTOR



LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE  
TL: 19 cm  
WL: 9 cm  
Tip Size: 1 mm  
30° Reverse Down  
Straight  
07-CVM-04



LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE  
TL: 19 cm  
WL: 9 cm  
Tip Size: 2 mm  
30° Reverse Down  
Straight  
07-CVM-05



LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE  
TL: 19 cm  
WL: 9 cm  
Tip Size: 3 mm  
30° Reverse Down  
Straight  
07-CVM-06



LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE  
TL: 19 cm  
WL: 9 cm  
Tip Size: 4 mm  
30° Reverse Down  
Straight  
07-CVM-07

Shown above at 100% scale  
Details shown larger

SI:3-2

DISSECTOR



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 1 mm  
30° Up  
Straight  
**07-CVM-08**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 2 mm  
30° Up  
Straight  
**07-CVM-09**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 3 mm  
30° Up  
Straight  
**07-CVM-10**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 4 mm  
30° Up  
Straight  
**07-CVM-11**

*Shown above at 100% scale  
Details shown larger*

SI-3-3

DISSECTOR

SI:3-4



LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE  
TL: 19 cm  
WL: 9 cm  
Tip Size: 1 mm  
45° Down  
Straight  
07-CVM-12



LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE  
TL: 19 cm  
WL: 9 cm  
Tip Size: 2 mm  
45° Down  
Straight  
07-CVM-13



LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE  
TL: 19 cm  
WL: 9 cm  
Tip Size: 3 mm  
45° Down  
Straight  
07-CVM-14



LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE  
TL: 19 cm  
WL: 9 cm  
Tip Size: 4 mm  
45° Down  
Straight  
07-CVM-15

Shown above at 100% scale  
Details shown larger

DISSECTOR



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 1 mm  
45° Up  
Straight  
**07-CVM-16**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 2 mm  
45° Up  
Straight  
**07-CVM-17**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 3 mm  
45° Up  
Straight  
**07-CVM-18**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 4 mm  
45° Up  
Straight  
**07-CVM-19**

*Shown above at 100% scale  
Details shown larger*

SI-3-5

DISSECTOR

SI:3-6



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 1 mm  
1 mm neck  
30° Reverse Down  
Straight  
**07-CVM-20**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 2 mm  
1 mm neck  
30° Reverse Down  
Straight  
**07-CVM-21**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 3 mm  
1 mm neck  
30° Reverse Down  
Straight  
**07-CVM-22**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 4 mm  
1 mm neck  
30° Reverse Down  
Straight  
**07-CVM-23**

*Shown above at 100% scale  
Details shown larger*



DISSECTOR



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 1 mm  
1 mm neck  
30° Up  
Straight  
**07-CVM-24**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 2 mm  
1 mm neck  
30° Up  
Straight  
**07-CVM-25**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 3 mm  
1 mm neck  
30° Up  
Straight  
**07-CVM-26**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 4 mm  
1 mm neck  
30° Up  
Straight  
**07-CVM-27**

*Shown above at 100% scale  
Details shown larger*

SI-3-7

DISSECTOR



SI:3-8



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 1 mm  
Straight  
**07-CVM-28**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 2 mm  
Straight  
**07-CVM-29**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 3 mm  
Straight  
**07-CVM-30**



**LAWTON DEEP  
CAVERNOUS  
MALFORMATION  
ROUND KNIFE**  
TL: 19 cm  
WL: 9 cm  
Tip Size: 4 mm  
Straight  
**07-CVM-31**

*Shown above at 100% scale  
Details shown larger*

## ENDONASAL DISSECTOR



**KELLY ENDOSCOPIC  
ENDONASAL  
DISSECTOR  
RHOTON #3 UP**  
TL: 22 cm  
WL: 12 cm  
Tip Size: 3 mm  
45° Up  
Bayonet  
**02-034-09MA**



**KELLY ENDOSCOPIC  
ENDONASAL  
DISSECTOR  
RHOTON #3 DOWN**  
TL: 22 cm  
WL: 12 cm  
Tip Size: 3 mm  
45° Down  
Bayonet  
**02-034-10MA**

*Shown above at 100% scale  
Details shown larger*

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

Part numbers may change without notice.  
For US commercial use only.  
Refer to the Instruction for Use before using.  
The colors in the photos may differ from the actual items.

[www.mizuho.com](http://www.mizuho.com)



©Mizuho America, ABPM 005, Rev A, Effective Date 09/28/2018

MUCOSAL DISSECTOR

SI:3-10



**KELLY ENDONASAL  
LEFT MUCOSAL  
DISSECTOR, SPEAR-LIKE**  
TL: 24 cm  
WL: 12.5 cm  
Spear-Like Shape  
Left  
Bayonet  
02-034-10



**KELLY ENDONASAL  
RIGHT MUCOSAL  
DISSECTOR, SPEAR-LIKE**  
TL: 24 cm  
WL: 12.5 cm  
Spear-Like Shape  
Right  
Bayonet  
02-034-11

*Shown above at 100% scale  
Details shown larger*

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

[www.mizuho.com](http://www.mizuho.com)



Part numbers may change without notice.  
For US commercial use only.  
Refer to the Instruction for Use before using.  
The colors in the photos may differ from the actual items.

## SKULL BASE DISSECTOR



### SKULL BASE DISSECTOR

A-1 5 mm Dissector  
TL: 19 cm  
WL: 9 cm  
Tip Size: 5 mm  
Sharp  
Slight Curve  
Round  
07-820-01



### SKULL BASE DISSECTOR

A-2 3.5 mm Dissector  
TL: 19 cm  
WL: 9 cm  
Tip Size: 3.5 mm  
Sharp  
Slight Curve  
Round  
07-820-02



### SKULL BASE DISSECTOR

A-3 2.5 mm Dissector  
TL: 19 cm  
WL: 9 cm  
Tip Size: 2.5 mm  
Sharp  
Slight Curve  
Round  
07-820-03



SI:3-12

*Shown above at 100% scale  
Details shown larger*

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

[www.mizuho.com](http://www.mizuho.com)



Part numbers may change without notice.  
For US commercial use only.  
Refer to the Instruction for Use before using.  
The colors in the photos may differ from the actual items.

©Mizuho America, ABPM 005, Rev A, Effective Date 09/28/2018

## SKULL BASE DISSECTOR



### SKULL BASE DISSECTOR

B-1 3.5 mm Dissector  
 TL: 19 cm  
 WL: 9 cm  
 Tip Size: 3.5 mm  
 Sharp  
 Slight Curve  
**07-820-04**



### SKULL BASE DISSECTOR

B-2 2.7 mm Dissector  
 TL: 19 cm  
 WL: 9 cm  
 Tip Size: 2.7 mm  
 Sharp  
 Slight Curve  
 Round  
**07-820-05**



### SKULL BASE DISSECTOR

B-3 2.5 mm Dissector  
 TL: 19 cm  
 WL: 9 cm  
 Tip Size: 2.5 mm  
 Sharp  
 Slight Curve  
**07-820-06**

SI:3-13

*Shown above at 100% scale  
 Details shown larger*

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

[www.mizuho.com](http://www.mizuho.com)



Part numbers may change without notice.  
 For US commercial use only.  
 Refer to the Instruction for Use before using.  
 The colors in the photos may differ from the actual items.

©Mizuho America, ABPM 005, Rev A, Effective Date 09/28/2018

## SKULL BASE DISSECTOR

---

SI:3-14



### SKULL BASE DISSECTOR

C-1 1.8 mm Dissector

TL: 19 cm

WL: 9 cm

Tip Size 1.8 mm

Sharp

Slight Curve

Round

07-820-07



### SKULL BASE DISSECTOR

C-2 1.8 mm Dissector

TL: 19 cm

WL: 9 cm

Tip Size 1.8 mm

Sharp

Moderate Curve

Round

07-820-08

*Shown above at 100% Scale*

*Details shown larger*

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

[www.mizuho.com](http://www.mizuho.com)



Part numbers may change without notice.  
For US commercial use only.  
Refer to the Instruction for Use before using.  
The colors in the photos may differ from the actual items.

©Mizuho America, ABPM 005, Rev A, Effective Date 09/28/2018

## SKULL BASE DISSECTOR



### SKULL BASE DISSECTOR

D-1 3.0 mm Dissector  
 TL: 19 cm  
 WL: 9 cm  
 Tip Size: 3.0 mm  
 Dull  
 Slight Curve  
 07-820-10



### SKULL BASE DISSECTOR

D-2 2.2 mm Dissector  
 TL: 19 cm  
 WL: 9 cm  
 Tip Size: 2.2 mm  
 Dull  
 Moderate Curve  
 Round  
 07-820-11



### SKULL BASE DISSECTOR

D-3 1.8 mm Dissector  
 TL: 19 cm  
 WL: 9 cm  
 Tip Size: 1.8 mm  
 Dull  
 Moderate Curve  
 Round  
 07-820-12

SI:3-15

*Shown above at 100% Scale  
 Details shown larger*

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

[www.mizuho.com](http://www.mizuho.com)



Part numbers may change without notice.  
 For US commercial use only.  
 Refer to the Instruction for Use before using.  
 The colors in the photos may differ from the actual items.

©Mizuho America, ABPM 005, Rev A, Effective Date 09/28/2018



# DISSECTORS

## SPATULA DISSECTOR



**KELLY ENDONASAL  
SPATULA DISSECTOR**

TL: 23.5 cm  
WL: 12.5 cm  
Tip Size: 2.4 mm  
Curved Up  
Bayonet  
02-034-12



**KELLY ENDONASAL  
SPATULA DISSECTOR**

TL: 23.5 cm  
WL: 12.5 cm  
Tip Size: 2.4 mm  
Curved Down  
Bayonet  
02-034-13

SI:3-16

*Shown above at 100% scale  
Details shown larger*

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

[www.mizuho.com](http://www.mizuho.com)



Part numbers may change without notice.  
For US commercial use only.  
Refer to the Instruction for Use before using.  
The colors in the photos may differ from the actual items.

©Mizuho America, ABPM 005, Rev A, Effective Date 09/28/2018



PANAMEDICAL SISTEMAS LTDA.  
Rua Borges Lagoa, 423  
Vila Clementino CEP 04038-030  
São Paulo/SP - Brasil  
Telefone 11 2344-1900  
[www.panamedical.com.br](http://www.panamedical.com.br)

IMPORTADOR E  
DISTRIBUIDOR  
EXCLUSIVO

**NSK**

 **MIZUHO**

 **FEATHER**®

 **TOITU**

**Mityvac**®



**CENTRO CIRÚRGICO E NEUROCIRURGIA**

[comercial@panamedical.com.br](mailto:comercial@panamedical.com.br)

11 2344-1903 | 11 94997-0139

**GINECO-OBSTETRÍCIA E MOBILIDADE**

[comercial.e@panamedical.com.br](mailto:comercial.e@panamedical.com.br)

11 2344-1905 | 11 98129-0112

**ASSISTÊNCIA TÉCNICA**

[atec@panamedical.com.br](mailto:atec@panamedical.com.br)

11 2344-1908 | 11 96437-0202