



PANAMEDICAL



LINHA DE PRODUTOS

DÚVIDAS E ORÇAMENTOS

comercial@panamedical.com.br | 11 23441900

Sugita AVM

MicroClips and Applier



5mm

4mm

3mm

2mm

Sugita AVM Microclips are made of cobalt-chromium-molybdenum-based alloy having superior properties, specified under ASTM F1058 G2(1S05832-7).

Being made of a non-ferromagnetic alloy, they are safe for MRI currently in clinical use.

To prevent damage and to facilitate handling, the clip is individually seated in a heat-resistant silicone sponge and packaged in a metal container, which is fully autoclavable.

Clips are made to have a holding force within the predetermined range of 50-70g.

- 07-932-02 AVM Clip, 2mm**
- 07-932-03 AVM Clip, 3mm**
- 07-932-04 AVM Clip, 4mm**
- 07-932-05 AVM Clip, 5mm**
- 07-935-10 AVM Clip Applier**
- 07-932-00 AVM Clip Case**



Clip Packaging



Mizuho America, Inc.
30057 Ahern Avenue
Union City, CA 94587

P: 800.699.2547
F: 510.324.4545
E: mizuho@mizuho.com
mizuho.com

