



PANAMEDICAL



LINHA DE PRODUTOS

DÚVIDAS E ORÇAMENTOS

comercial@panamedical.com.br | 11 23441900

SI 8 RING CURETTE

For Ordering information contact your Sales Rep
or Mizuho America at 800-699-CLIP (2547)
www.MIZUHO.com

CP ANGLE RING CURETTE



CP-8 TUMOR RING CURETTE
 TL: 19.5 cm
 WL: 9 cm
 Tip Size: 1.2 mm
 60°
 Round
 07-821-08



CP-9 TUMOR RING CURETTE
 TL: 19.5 cm
 WL: 9 cm
 Tip Size: 1.55 mm
 45°
 Round
 07-821-09



CP-10 CUP CURETTE
 TL: 19.5 cm
 WL: 9 cm
 Tip Size: 1 x 1.5 mm
 45°
 Round
 07-821-10



CP-11 CUP CURETTE
 TL: 19.5 cm
 WL: 9 cm
 Tip Size: 1.5 x 2 mm
 45°
 Round
 07-821-11

Shown above at 100% scale
 Details shown larger

RING CURETTE

SI:8-2



**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**

TL: 23.5 cm
WL: 12.5 cm
Tip Size: 1.5 mm
45° Up
Bayonet
07-822-21



**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**

TL: 23.5 cm
WL: 12.5 cm
Tip Size: 3 mm
45° Up
Bayonet
07-822-23



**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**

TL: 23.5 cm
WL: 12.5 cm
Tip Size: 5 mm
45° Up
Bayonet
07-822-25



**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**

TL: 23 cm
WL: 12.5 cm
Tip Size: 7 mm
45° Up
Bayonet
07-822-27

Shown above at 100% scale
Details shown larger

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

Part numbers may change without notice.
For US commercial use only.
Refer to the Instruction for Use before using.
The colors in the photos may differ from the actual items.

www.mizuho.com



RING CURETTE



⊙ Actual Size

**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**

TL: 23 cm
WL: 12.5 cm
Tip Size: 1.5 mm
90° Up
Bayonet
07-822-31



⊙ Actual Size

**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**

TL: 23 cm
WL: 12.5 cm
Tip Size: 3 mm
90° Up
Bayonet
07-822-33



⊙ Actual Size

**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**

TL: 23 cm
WL: 12.5 cm
Tip Size: 5 mm
90° Up
Bayonet
07-822-35



⊙ Actual Size

**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**

TL: 23 cm
WL: 12.5 cm
Tip Size: 7 mm
90° Up
Bayonet
07-822-37

*Shown above at 100% scale
Details shown larger*

SI-8-3

RING CURETTE

SI:8-4



**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**
TL: 23 cm
WL: 12.5 cm
Tip Size: 1.5 mm
45° Down
Bayonet
07-822-41



**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**
TL: 23 cm
WL: 12.5 cm
Tip Size: 3 mm
45° Down
Bayonet
07-822-43



**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**
TL: 23 cm
WL: 12.5 cm
Tip Size: 5 mm
45° Down
Bayonet
07-822-45



**KELLY ENDOSCOPIC
ENDONASAL
RING CURETTE**
TL: 23 cm
WL: 12.5 cm
Tip Size: 7 mm
45° Down
Bayonet
07-822-47

Shown above at 100% scale
Details shown larger

RING CURETTE



KELLY ENDOSCOPIC ENDONASAL RING CURETTE

TL: 23 cm
WL: 12.5 cm
Tip Size: 7 mm
Left to Right Action
Bayonet
07-822-97

*Shown above at 100% scale
Details shown larger*

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

www.mizuho.com



Part numbers may change without notice.
For US commercial use only.
Refer to the Instruction for Use before using.
The colors in the photos may differ from the actual items.

©Mizuho America, ABPM 005, Rev A, Effective Date 09/28/2018

SI-8-5

RING CURETTE

RING CURETTE

SI:8-6



WL: 9 cm
WL: 11 cm

KELLY KEYHOLE RING CURETTE

Tip Size: 1.5 mm
45° Up
Bayonet

07-822-60MA
07-822-63MA



⊙ Actual Size

KELLY KEYHOLE RING CURETTE

Tip Size: 3 mm
45° Up
Bayonet

07-822-54MA
07-822-57MA



⊙ Actual Size

KELLY KEYHOLE RING CURETTE

Tip Size: 5 mm
45° Up
Bayonet

07-822-48MA
07-822-51MA



⊙ Actual Size

*Shown above at 100% scale
Details shown larger*

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

www.mizuho.com



Part numbers may change without notice.
For US commercial use only.
Refer to the Instruction for Use before using.
The colors in the photos may differ from the actual items.

©Mizuho America, ABPM 005, Rev A, Effective Date 09/28/2018

RING CURETTE



WL: 9 cm
WL: 11 cm

**KELLY KEYHOLE
RING CURETTE**

Tip Size: 1.5 mm
90° Up
Bayonet

07-822-61MA
07-822-64MA



**KELLY KEYHOLE
RING CURETTE**

Tip Size: 3 mm
90° Up
Bayonet

07-822-55MA
07-822-58MA



**KELLY KEYHOLE
RING CURETTE**

Tip Size: 5 mm
90° Up
Bayonet

07-822-49MA
07-822-52MA



*Shown above at 100% scale
Details shown larger*

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

www.mizuho.com



Part numbers may change without notice.
For US commercial use only.
Refer to the Instruction for Use before using.
The colors in the photos may differ from the actual items.

©Mizuho America, ABPM 005, Rev A, Effective Date 09/28/2018

SI-8-7

RING CURETTE

RING CURETTE

SI:8-8



WL: 9 cm
WL: 11 cm



© Actual Size

KELLY KEYHOLE RING CURETTE

Tip Size: 1.5 mm
45° Down
Bayonet

07-822-62MA
07-822-65MA



© Actual Size

KELLY KEYHOLE RING CURETTE

Tip Size: 3 mm
45° Down
Bayonet

07-822-56MA
07-822-59MA



© Actual Size

KELLY KEYHOLE RING CURETTE

Tip Size: 5 mm
45° Down
Bayonet

07-822-50MA
07-822-53MA

Shown above at 100% scale
Details shown larger

For Ordering information contact your Sales Rep or Mizuho America at 800-699-CLIP (2547)

www.mizuho.com



Part numbers may change without notice.
For US commercial use only.
Refer to the Instruction for Use before using.
The colors in the photos may differ from the actual items.

©Mizuho America, ABPM 005, Rev A, Effective Date 09/28/2018



PANAMEDICAL SISTEMAS LTDA.
Rua Borges Lagoa, 423
Vila Clementino CEP 04038-030
São Paulo/SP - Brasil
Telefone 11 2344-1900
www.panamedical.com.br

IMPORTADOR E
DISTRIBUIDOR
EXCLUSIVO

NSK

 **MIZUHO**

 **FEATHER**®

 **TOITU**

Mityvac®



CENTRO CIRÚRGICO E NEUROCIRURGIA

comercial@panamedical.com.br
11 2344-1903 | 11 94997-0139

GINECO-OBSTETRÍCIA E MOBILIDADE

comercial.e@panamedical.com.br
11 2344-1905 | 11 98129-0112

ASSISTÊNCIA TÉCNICA

atec@panamedical.com.br
11 2344-1908 | 11 96437-0202